



## V4c Clinical Outcomes


**“The most impressive outcome was the intraocular pressure post-operatively. It was identical and normal in both sets of eyes. We have continued to use the technology in our patients since the trials and they are very happy.”**

Prof. Kimiya Shimizu, Japan


### Superb Quality of Vision\*



### IOP Stability\*



### Increased Surgical Efficiency\*



\* Data on file with STAAR Surgical  
 \*\* Igarashi A, Kamiya K, Shimizu K, Komatsu M. Visual Performance after Implantable Collamer Lens Implantation and Wavefront-Guided Laser In Situ Keratomileusis for High Myopia. Am J Ophthalmol. 2009; 148 (1):164-170

## Product Specification

### SPHERICAL LENSES

Lens Powers	Optical Zone Diameter (in mm)	Effective OZ at Corneal Plane (in mm)
-0.5 to -9.0	5.8	7.30
-9.5 to -10.0	5.5	6.93
-10.5 to -12.5	5.3	6.62
-13.0 to -18.0	4.9	6.17



0.25D increments available from -0.5D to -3.0D      Lens lengths: 12.1 mm / 12.6 mm / 13.2 mm / 13.7 mm

### TORIC LENSES

Lens Powers	Cylinder	Optical Zone (in mm)	Effective OZ at Corneal Plane (in mm)
-0.5 to -9.0	+0.5 to +6.0	5.8	7.30
-9.5 to -10.0	+0.5 to 6.0	5.5	6.93
-10.5 to -12.5	+0.5 to 6.0	5.3	6.62
-13.0 to -18.0	+0.5 to 6.0	4.9	6.17

Lens lengths: 12.1 mm / 12.6 mm / 13.2 mm / 13.7 mm


### ICL™ / Toric ICL™ Calculation interface



### Customized Implantation Diagrams for Each Patient

The ICL™ Power Calculation Tool provides surgeons the ability to calculate the proper diopter, cylinder and lens size for given pre-operative measurements. STAAR targets to provide Toric ICLs with a customized cylinder axis allowing for easy horizontal lens alignment. Some lenses may be provided with minor axis deviations for which the implantation orientation diagram will assist in accurate lens alignment.

Distributed By




STAAR® Surgical Head Office · 1911 Walker Ave. · Monrovia, CA 91016 · +1 (800) 352-7842  
 STAAR® Surgical International Office · Hauptstrasse 104 · CH 2560 Nidau/Switzerland · +41 32 332 8888  
 www.staar.com

Made in Switzerland



# CLEARLY MORE SIMPLE



## The New Visian ICL® with CentraFLOW™ Technology

## The new Visian ICL<sup>®</sup> V4c with CentraFLOW<sup>™</sup> Technology clearly simplifies the ICL surgical procedure

- Eliminates the need for PIs ; increasing the efficiency for both surgeon and patient
- Restores a more natural aqueous flow
- Facilitates OVD removal

**“This revolutionary CentraFLOW<sup>™</sup> Technology will be a game changer for implant-based refractive surgery”**

Dr. Erik Mertens, Belgium


**“With the V4c, there is no longer a fear of PI failure. You don't have to worry about pupillary block. And for the patient, having no PI means quicker surgery and faster visual recovery.”**

Dr. Daniel Elies, Spain

### CentraFLOW<sup>™</sup> Technology

- Efficacious
- IOP Stability
- Superb Quality of Vision
- Increased Surgical Efficiency
- Enhanced Convenience and Comfort for the Patient

## Over a quarter of a million ICLs implanted worldwide



KS-Aquaport<sup>™</sup>


## The KS-Aquaport<sup>™</sup>

**“The KS-AP technology was aberration free and all visual results were excellent and identical in both the KS-AP and the conventional ICL groups.”**

Prof. Kimiya Shimizu, Japan


### How will an optic port affect quality of vision

The MTF values with and without a hole are nearly identical



### Optimizing the size of the port

A significant drop in optical quality occurs only with a hole >500μ



### Simulated Cross-Sectional Image of the KS-Aquaport



**Superior edge finish supports high quality of vision**

The highly polished edges of the port result in rounded smooth surfaces having virtually no impact on quality of vision.

\* Data on file with STAAR Surgical