

# IOLS



Reference	Model	Type	Filter	Brand name	Optical diameter	Total length	Diopters range	Edge	Water content	Sterilization method	Expiration	A Constant / ACD
P651300	POSTERIOR CHAMBER	Polymethyl methacrylate (PMMA)	UV	AJL posterior chamber IOL	6.5mm	13mm	-10 to +40 in steps of 1D / +12 to +25 in steps of 0.5D	Square	N/A	ETO	4.5 years	A Constant (U/S): 118 / ACD:4.96

# IOLS



Reference	Model	Type	Filter	Brand name	Optical diameter	Total length	Diopters range	Edge	Water Content	Sterilization method	Expiration	A Constant / ACD
A601250	ANTERIOR CHAMBER	Polymethyl methacrylate (PMMA)	UV	AJL anterior chamber IOL	6.0mm	12.5mm	+6 to +30D in steps of 1D	Square	N/A	ETO	4.5 years	A Constant (U/S):115,5 / ACD: 3,51