



The future of vision



INTACS[®] Intrastromal Corneal Implants

<i>Model</i>	Intacs (optical area of 7 mm) and Intacs SK (optical area of 6 mm).
<i>Brand name</i>	INTACS [®] and INTACS [®] SK
<i>Definition</i>	Intrastromal Corneal Implants
<i>Indications</i>	<p>Myopia: Intacs corneal implants are indicated for the reduction or correction of mild to moderate myopia (from -0.50 D to -5.00 D) in patients over 21 years, and astigmatism with diopters less than or equal to +1 D</p> <p>Keratoconus: Intacs corneal implants are indicated for the treatment of keratoconus, to fully or partially correct myopia or astigmatism in patients that need to restore their visual function and for whom a satisfactory correction of vision with contact lenses or with glasses is not possible, and for patients who wish to avoid a possible corneal transplant.</p>
<i>Product description</i>	Intacs[®] and Intacs SK[®] intracorneal implants are rounded segments designed to implant in the corneal stroma with a small radial incision, in order to reduce halos and other undesired visual effects.
<i>Material</i>	Polymethyl Methacrylate (PMMA)
<i>Sterilization method</i>	Ethylene Oxide (ETO)
<i>Supply</i>	Sterile, single or double units, depending on the reference requested (<i>see table</i>).
<i>Shelf-life</i>	5 years

CE 0086



REV: 01



INTACS CORNEAL IMPLANT	
<i>PRODUCT</i>	<i>SINGLE SEGMENT (1 UNIT)</i>
DESCRIPTION	Intacs® (7.0 mm)
Corneal Implant 0.250 mm-210°	12400-250-210
<i>PRODUCT</i>	<i>SYMMETRICAL TREATMENT (2 UNITS)</i>
DESCRIPTION	Intacs® (7.0 mm)
Corneal Implant 0.210 mm-150°	11500-210-150
Corneal Implant 0.250 mm-150°	11500-250-150
Corneal Implant 0.300 mm-150°	11500-300-150
Corneal Implant 0.350 mm-150°	11500-350-150
Corneal Implant 0.400 mm-150°	11500-400-150
Corneal Implant 0.450 mm-150°	11500-450-150
<i>PRODUCT</i>	<i>ASYMMETRICAL TREATMENT (2 UNITS)</i>
DESCRIPTION	Intacs® (7.0 mm)
Corneal Implant 0.210/0.300-150°	11200-210/300-150
Corneal Implant 0.210/0.350-150°	11200-210/350-150
Corneal Implant 0.210/0.400-150°	11200-210/400-150
Corneal Implant 0.210/0.450-150°	11200-210/450-150

*OTHER THICKNESSES AVAILABLE UPON REQUEST



INTACS SK CORNEAL IMPLANT		
<i>PRODUCT</i>	<i>SINGLE SEGMENT (1 UNIT)</i>	<i>SYMMETRICAL TREATMENT (2 UNITS)</i>
DESCRIPTION	Intacs® SK (6.0 mm)	Intacs® SK (6.0 mm)
Corneal Implant 0.250 mm-90°	12300-250-90	
Corneal Implant 0.300 mm-90°	12300-300-90	
Corneal Implant 0.350 mm-90°	12300-350-90	
Corneal Implant 0.400 mm-90°	12300-400-90	
Corneal Implant 0.450 mm-90°	12300-450-90	
Corneal Implant 0.300 mm-130°	12300-300-130	
Corneal Implant 0.350 mm-130°	12300-350-130	
Corneal Implant 0.400 mm-130°	12300-400-130	
Corneal Implant 0.450 mm-130°	12300-450-130	
Corneal Implant 0.210 mm-150°	12300-210-150	12000-210-150
Corneal Implant 0.250 mm-150°	12300-250-150	12000-250-150
Corneal Implant 0.300 mm-150°	12300-300-150	12000-300-150
Corneal Implant 0.350 mm-150°	12300-350-150	12000-350-150
Corneal Implant 0.400 mm-150°	12300-400-150	12000-400-150
Corneal Implant 0.450 mm-150°	12300-450-150	12000-450-150
Corneal Implant 0.500 mm-150°	12300-500-150	12000-500-150

*OTHER THICKNESSES AVAILABLE UPON REQUEST