



OPTIM[®]33TB Wipes Quick, Clean and Green!

- ✔ - Superior quality design
- ✔ - Non woven
- ✔ - Non abrasive
- ✔ - Economical
- ✔ - Broad spectrum
- ✔ - Tuberculocidal, Virucidal, Bactericidal




Superior quality wipes for surface cleaning & disinfection.

50 YEARS OF
INNOVATION
1957-2007

SciCan
A HIGHER STANDARD

OPTIM 33TB Wipes – Superior quality design



Unique, non woven bond pattern with uniform pore size ensures that the liquid is present on the substrate rather than being absorbed. Result:

- The liquid is deposited onto the surface instead of staying inside the wipe
- Optim wipes remain extremely wet and do not dry out
- Durable and strong with no lint residue
- Capable of covering a large area with minimal wipe use
- Economical consumption with less waste

OPTIM 33TB Wipes – Quick

- **30 Second sanitizer** - 99.999% broad spectrum kill (5-log reduction) of vegetative bacteria
- **1 Minute** Tuberculocidal, Bactericidal, Virucidal (broad spectrum - including hard to kill non-enveloped viruses)
- Other Microbial Data on file: 3 minute Fungicide†

OPTIM 33TB Wipes – 2 formats available



- **Ready to use presaturated wipes**

SCI33W12-ALOE Optim 33 TB RTU Wipes with Aloe Scent 6 x 7"
(1 case: 12 x 160 sheets)



- **Dry wipes in two sizes allow for reuse of Optim wipes canister in combination with Optim liquid. Ideal for multi-operatory practices providing cost savings and waste reduction.**

Dry Wipes:

SCI33R18 Optim 33 TB Dry Wipes Refill 6 x 7"
(1 case: 18 x 160 sheets)

SCI33R18XL Optim 33 TB Large Dry Wipes Refill 10 x 10"
(1 case: 18 x 60 sheets)

Liquid:

SCI3312X1-ALOE Optim 33 TB with Aloe Scent, 1L (1 case: 12 x 1L)

SCI334X4 Optim 33 TB unscented, 4L (1 case: 4 x 4L)

SCI334X4-ALOE Optim 33 TB with Aloe Scent, 4L (1 case: 4 x 4L)

SciCan

A HIGHER STANDARD

SciCan Ltd., 1440 Don Mills Road, Toronto, Ontario M3B 3P9
Phone (416) 445 1600 Fax (416) 445-2727 Toll free 1-800-667-7733

www.scican.com

