


# Introduction to the AcuTarget™ Diagnostic System

**AcuFocus™**

# AcuTarget™ System Overview

## Diagnostic Unit




# Benefits at a Glance


- Simplify & speed up workflow
- Reliable pre-op & post-op diagnostics
- Reproducible precision
- Pre & post-op documentation

# Pre-Operative Functionality

- Identify:
  - 1<sup>st</sup> Purkinje Reflex
  - Pupil boundary and centroid
- Inlay centration location & planning




# Pre-Op Captured Information




# Post-Op Functionality

- Identify the inlay position
  - Relative to 1st Purkinje
  - Relative to Pupil center
  - Pupil boundary
  - Limbus



# Post-Op Captured Information



# Summary

- Precisely and reliably centers the inlay
- Improves patient outcomes
- Reduces the need for recentration
- Reduces exam/chair time





Thank You