



OCULAR CONFORMERS



Ocular Conformers are designed to be used after surgery to prevent closure or adhesions during the healing process. Conformers are small, acrylic-cup shaped devices whose inner surfaces are shaped to approximate the curvature of the orbit.

MEDPOR Ocular Conformers are supplied sterile.



MEDPOR Ocular Conformers

CAT#	DESCRIPTION	A	B	C
9778	Vented - Small	19mm	x 22mm	x 4.4mm
9779	Vented - Med	21mm	x 24mm	x 5.3mm
9780	Vented - Large	22mm	x 25mm	x 5.6mm