



EYE SHIELD

The future of vision



AJL Eye Shield

<i>Brand name</i>	AJL EYE SHIELD		
<i>Reference</i>	PROCUL-S-AJL	PROCUL-M-AJL	PROCUL-L-AJL
<i>Size</i>	Small: 22 x 19 x 1.7 mm	Medium: 24 x 21 x 1.7 mm	Large: 25 x 22 x 1.7 mm
<i>Product description</i>	AJL EYE SHIELD is a medical device specially designed to protect the eyeball.		
<i>Material / Composition</i>	Acrylic		
<i>Indications</i>	AJL EYE SHIELD is designed to protect the eyeball during eyelid surgery or in DCR (Dacriocystorriostomy).		
<i>Sterilization method</i>	Ethylene Oxide (ETO)		
<i>Supply</i>	AJL EYE SHIELD is presented sterile in an easy opening holder. It is provided together with a suction teat that helps positioning and removal of the eye shield.		
<i>Expiration</i>	4.5 years		

CE 0434

REV:00