







LIOCAN

AJL manufactures VET intraocular lenses based on the data obtained in eye specific study of different races, thus fixing the most optimized design for each one of them

VET	IOL MODEL (mm) 	DIOPTERS
 <p>Mastin Bóxer</p>	Ref: TXY651525 OPT: 6.5 OT: 15,25 Aspheric IOL	41
 <p>Perdigueros Labradores Lebreles</p>	Ref: TXY651425 OPT: 6.5 OT: 14,25 Aspheric IOL	41
 <p>Cocker Spaniel Foxterrier Beagle</p>	Ref: TXY651325 OPT: 6.5 OT: 13,25 Aspheric IOL	41
 <p>Teckel WHW Terrier Eichón Maltés</p>	Ref: TXY651225 OPT: 6.5 OT: 12,25 Aspheric IOL	41
 <p>Yorkshire Terrier Carlino Fomerania</p>	Ref: TXY651125 OPT: 6.5 OT: 11,25 Aspheric IOL	41

An injection system is included. Medicell Naviject injector 2.9mm-1P

www.ajlsa.com