



SOLEKO



ALL IN ONE SOLUTION

HYDRATING DROPS

YAL PRODUCT LINE

ALL IN ONE SOLUTION

YALFRESH

FEATURES

YALFRESH All in One Solution is characterized by its formulation based on **Hyaluronic acid**, a natural polysaccharide capable of increasing the wettability of the lenses, prolonging the comfort of contact lens wearers.

Hyaluronic acid has excellent hygroscopic properties. In contactology, this extraordinary capacity for water retention is used to improve the hydration of the precorneal area. Another important feature of **Hyaluronic acid** is mucus-adhesiveness, which allows it to form a long-lasting coating and provides a stable protection to the corneal surface. Finally, it has viscoelastic properties that lubricate the ocular surface, reducing surface friction during eye movements and eyelid blinking.

YALFRESH All in One Solution is indicated for the disinfection, storage, cleansing, rinsing, hydration and prevention of protein deposits of all types of soft and color contact lenses.

INGREDIENTS

Ethylene diamine tetra acetic acid 0.01%, PHMB 0.00015%, Hydroxyethylcellulose, Hyaluronic acid macromolecular, Natrium chloride, Dibasic sodium phosphate, Monobasic sodium phosphate, Sodium hydroxide anydrus, Purified water.



ALL IN ONE SOLUTION

YALFRESH

BENEFITS

Effectiveness in disinfection and film-forming, emollient and nourishing activity.
It respects and protects the corneal epithelium.

Efficiency in cleaning contact lenses.
Extended comfort during contact lens wear.

Formulation containing hyaluronic acid.
It prevents dehydration of contact lenses.



Compatibile con i tessuti ed i liquidi oculari
Compatible with ocular tissues/fluids



Deterge
Cleans



Idrata
Hydrates



Umetta
Lubricates



Previene i depositi proteici
Prevents protein deposits



Disinfetta
Disinfects



Risciacqua
Rinse



Indicato per lenti morbide
Suitable for soft contact lenses



Acido Ialuronico
Hyaluronic Acid

I-FRESHYAL

FEATURES

The humectant solution **I-FRESH YAL** is characterized by its formulation containing **Hyaluronic acid**.

The lubricating and viscosifying properties of Hyaluronic acid reduce the friction generated by the rubbing of the eyelid-contact lens-cornea, prevent and soothe the typical disorders related to the tear film instability and alleviate ocular irritations caused by mechanical stress due to a foreign body or after the use of any type of contact lens.

I-FRESH YAL solution restores the hydration of soft lenses and the wettability of rigid gas permeable ones. It is indicated as a lubricant solution during contact lens application, in case of tired and poor tearing eyes, or when sensations of foreign bodies are felt into the eyes.



INGREDIENTS

Ethylene diamine tetra acetic acid 0.023%, PHMB 0.00005%, Hydroxyethylcellulose, Hyaluronic acid macromolecular 0.2%, Natrium chloride, Dibasic sodium phosphate, Monobasic sodium phosphate, Sodium hydroxide anydrus, Purified water.

HYDRATING DROPS

I-FRESHYAL

BENEFITS

Humectant action on contact lenses.

Extended comfort during contact lens wear.

Hydrating and refreshing action.

It relieves burning symptoms also due to poor blinking.

It stabilize the tear film.

It contributes to the maintenance of ocular homeostasis.



Compatibile con i tessuti ed i liquidi oculari
Compatible with ocular tissues/fluids



Idrata
Hydrates



Umetta
Lubricates



Indicato per lenti morbide
Suitable for soft contact lenses



Silicone Hydrogel
Silicone hydrogel



Indicato per lenti RGP
Suitable for RGP contact lenses

