



Orbital Implant Oculfit

OCULFIT

Reference	POI1400	POI1600	POI1800	POI1900	POI2000	POI2100	POI2200	POI2300
Size	14mm Diameter Sphere	16mm Diameter Sphere	18mm Diameter Sphere	19mm Diameter Sphere	20mm Diameter Sphere	21mm Diameter Sphere	22mm Diameter Sphere	23mm Diameter Sphere
Model	Section "Reference"							
Brand name	ORBITAL IMPLANT "OCULFIT"							
Definition	Porous polyethylene spheres. Orbital implants.							
Indications	The orbital implants OCULFIT are spherical implants which are designed to be implanted into resulting cavity of eviscerated or enucleated eyeball							
Product description	<p>Due to the structure of interconnected opened porous, this product allows the growth of the tissue which is inside it.</p> <p>OCULFIT orbital implants from AJL have smooth porous anterior surface which helps to minimize the expose of the implants for a long time of period.</p> <p>Biocompatibility studies in vitro and in vivo have demonstrated that OCULFIT implants are free from cytotoxic effects.</p>							
Size	Section "Reference"							
Material & Composition	High density Polyethylene. According to different studies and experiences, this material is well tolerated by the tissues and it has a good stability.							
Presentation	OCULFIT orbital implants are sterile and non-pyrogenic. The product is packed inside a Tyvek bag that maintained its sterility during all implant's useful life. There are sterilized for Ethylene Oxide.							
Sterilization method	Ethylene Oxide (ETO)							
Expiration	5 years							

